


Products


Product	CAT#	Purity	Source	Application
SAA1 Ag	FAP-BE004	>95%	E.coli	Calibrator/ Quality Control

Product	CAT#	Clone	Source	Pair/ Application
SAA Ab	FAB-B004-C1	C1	CHO	C1-C2 (CLIA) C2-C1 (CLIA) C3-C4 (TRFIA) 2C5 (LETIA)
	FAB-B004-C2	C2	CHO	
	FAB-B004-C3	C3	CHO	
	FAB-B004-C4	C4	CHO	
	FAB-B004-2C5	2C5	Mouse	

Clinical Performance



Clinical Performance



Detection of Siemens assigned clinical serum
(N=11, Concentration range: 0.2-110 μ g/ml)
Sensitivity: 7ng/ml,
Sample compliance rate: $R^2 > 95\%$.

Contact us:
3rd F, MTP, Shuitu , Liangjiang New District,
Chongqing China 400714
Tel: +86-23-68616605
Email: info@fantibody.com
website: www.fantibody.com

Quotes:
sales@fantibody.com
Shipment status:
order@fantibody.com
Technical support:
support@fantibody.com

