


Products

Product	CAT#	Purity	Source	Application
AMH Ag	FAP-PC003	>80%	293 cells	Calibrator or Quality Control
	FAP-PC003N	>95%	293 cells	Calibrator or Quality Control

Product	CAT#	Clone	Recommended use	Pair/ Application
AMH Ab	FAB-P003	3D3	Capture/ Detection	1B9 - 3D3
		1B9	Capture/ Detection	
		6B3	Capture	6B3 - 1B9
		2C10	Detection	3D3 - 2C10
		3G8	Capture	3G8-3G10
		3G10	Detection	


Clinical Performance




AMH linear range

Sensitivity : below 0.05ng/mL
Linear range: 0.05-20ng/mL

Clinical Performance



AMH clinical relevance $R^2 = 0.9648$




AMH linear range

Sensitivity : below 0.05ng/mL

Linear range: 0.05-20ng/mL

Clinical Performance



AMH clinical relevance $R^2 = 0.9768$

Contact us:
3rd F, MTP, Shuitu , Liangjiang New District,
Chongqing China 400714
Tel: +86-23-68616605
Email: info@fantibody.com
website: www.fantibody.com

Quotes:
sales@fantibody.com
Shipment status:
order@fantibody.com
Technical support:
support@fantibody.com

