

# Profile of PCT

## Products

Product	CAT#	Purity	Source	Application
PCT Ag	FAP-BE002	>95%	E.coli	Calibrator or Quality Control
PCT Ag	FAP-BC002	>95%	CHO cells	Calibrator or Quality Control


Product	CAT#	Clone	Source	Capture/ Detection	Pair/ Application
PCT Ab	FAB-B002	10D6	Mouse	Capture	10D6-4A1
		4A1	Mouse	Detection	

## Clinical Performance


### PCT antigen stability

Temperature	time	0.1 ng/ml	1 ng/ml	10 ng/ml	60 ng/ml
-20℃	1周	1.428	2.027	2.501	2.508
37℃	1周	0.97	1.671	2.122	2.248
PCT antigen , 10D6 (Capture) - 4A1 (Detection) sandwich					

# Clinical Performance



PCT antibody activity, direct ELISA  
PCT antigen 2ug/well capture



10D6-4A1 linear range

Contact us:  
3rd F, MTP, Shuitu , Liangjiang New District,  
Chongqing China 400714  
Tel: +86-23-68616605  
Email: info@fantibody.com  
website: www.fantibody.com

Quotes:  
sales@fantibody.com  
Shipment status:  
order@fantibody.com  
Technical support:  
support@fantibody.com

