

Profile of HBP


Products

Product	CAT#	Purity	Source	Application
HBP Ag	FAP-BC005	> 95%	HEK 293F cells	Calibrator/ Quality Control

Product	CAT#	Clone	Source	Capture/ Detection	Pair/ Application
HBP Ab	FAB-B005	10G11	Mouse	Detection	19A8-10G11 13D8-12F5 13D8-10G11
		12F5	Mouse	Detection	
		13D8	Mouse	Capture	
		19A8	Mouse	Detection	

Clinical Performance

HBP Antigen purity





Antigen activity assay: capture 1:1000, HRP 1:9000

antigen dilution	1:5k	1:10k	1:20k	1:40k	1:80k	1:160k	1:320k	0
Signal value	1147288	1057806	1049805	772588	606099	432483	244132	32111


Antigen stability test: capture 1:2000, HRP 1:1000

antigen dilution	1:20k	1:40k	1:80k	1:160k	1:320k	1:640k	1:1280k	0
initial value	1614566	1249842	980973	472185	253346	175574	107801	16572
37°C 3days	1607133	1284640	880612	616539	333098	204184	128304	32840

Pairing information



Pairing information



LFIA: sensitivity: 0.5ng/ml,
Linear range: 0.6-300ng/ml

Contact us:
3rd F, MTP, Shuitu , Liangjiang New District,
Chongqing China 400714
Tel: +86-23-68616605
Email: info@fantibody.com
website: www.fantibody.com

Quotes:
sales@fantibody.com
Shipment status:
order@fantibody.com
Technical support:
support@fantibody.com

